


**Legacy 8 Portable Basketball Goal**  
Overhead Collapsed View (Anchor Hidden)



14'1" is the amount of clear space needed for a standard Kinematic Integrated Stabilizing System (KISS) Legacy Portable Basketball Goal. This is the space needed to position the goal behind the baseline while in the collapsed position.

Only 10'8" of space (10'11" to the rear of the anchor cover if applicable) is needed when the goal is in the playing position.

The Manual KISS version of Legacy should be ordered if you lack the additional 3'5" of clear space. This allows the goal to be raised into the playing position first, and then rolled into place behind the baseline.